



The Invisiport™ is a micro-invasive implantable port offering substantial benefits to patients and medical professionals by providing port convenience with PICC simplicity. The Invisiport™ is the only port that can be placed through a 10 mm dermatotomy.

INVISIPORT™

Specifications for the Invisiport™

The following table provides the specifications for the Invisiport™.

Catheter Material	polyurethane
Septum Material	titanium
Catheter Length	53 cm (20.87 inches) Catheter is trimmable to desired length
Catheter Outer Diameter	2.01 mm (6 Fr) (.094 inches)
Catheter Inner Diameter	1.30 mm (.0625 inches)
Power Injectable	300 psi/5 cc/sec

Accessory Items:

- One (1) Galt Microaccess Tear-away Introducer (6Fr x 50 cm/19.69 inches)
- One (1) Galt Microintroducer Needle
- One (1) Galt Microguidewire (Nitinol)
- One (1) 20 gauge, 1.5 inch, Non-Coring Access Needle

Product Number: 010966

For More Information

Stealth Therapeutics, Inc.
406 Science Drive, Suite 400
Madison, WI 53711
(877) 262-4946
www.stealththerapeutics.com
customerservice@stealththerapeutics.com

The Invisiport™ is indicated for any patient requiring repeated access of the vascular system for delivery of medications, nutritional supplementation, fluids, blood, blood products, or the sampling of blood, or for power injection of contrast media when used with a power-injectable Huber needle or infusion set.



Benefits of Invisiport™

Versus PICCs

- ▶ Improved quality of life
- ▶ More conducive to bathing and life activities
- ▶ Eliminate risk of dislodgement
- ▶ Potential for reduced infection¹

Versus Traditional Chest Ports

- ▶ Less anesthesia
- ▶ Shorter procedure time
- ▶ Smaller incision, smaller scar, less visible scar



1 - "Long-Term Central Venous Access," Stephen P. Povoski, M.D., November 2, 2011